

CONTROLS

Test to check sterilization and to detect pyrogenic substances, control to check flow and seal. Chemical control of components in accordance with the European Pharmacopoeia (Current edition).

STERILIZATION

Product sterilized with ethylene oxide (E.O.). Validity: 5 years from sterilization date.

FIDUUCI SIEI

STORAGE

Normal storage procedures.

DISPOSAL

Please follow the current laws on hospital waste disposal.

00	03/06/2020	Emission	L. Rossi	M. Mancarti	G. Galli
REV	DATE	DESCRIPTION	ISSUED BY	APPROVED BY	QA APPROVED
01	09/02/2021	Changed mould	L. Rossi	M. Mancarti	G. Galli

C € 0051

Via dell'industria 3 - 23030 Tovo di S. Agata (SO) Italy - Tel +39 0342. 771049 Fax +39 0342. 771709- info@valmedsrl.com