


## TPN Multilayer Bag 3000ml female connector

**CODE:** MB-1430


**PCS/BOX:** 50

**WEIGHT:** 100g

1. Renolit Solmed Medipak UVO 9002 400 micron
2. Internal handle (PS)
3. Black graduated scale
4. Coex Ø4.8x6.8mm tubes (PUR-EVA)
5. Coex Ø6X8mm tube (PUR-EVA)
6. Infusion port (ABS)
7. Safety cap for infusion port (LDPE)
8. Inviolable clamp (SBC)
9. Spike port (PVC DEHP free)
10. Female connector (ABS)
11. Male vented cap (ABS)
12. Male cap blister packaging (ABS)
13. Paper/PE envelope



**DEHP FREE**



### Intended use

Container for nutrition solutions prior and during administration to a patient.

### Controls

Test to check sterilization and to detect pyrogenic substances, control to check flow and seal.

Chemical control of components in accordance with the European Pharmacopoeia (Current edition).

### Sterilization

Product sterilized with ethylene oxide (E.O.). Validity: 5 years from sterilization date.

### Packing

Single package product in a Paper/PE envelope.

The product is supplied sterile and non pyrogenic.

### Storage

Normal storage procedures.

### Disposal

Please follow the current laws on hospital waste disposal.

**Miktell Ltd**  
 Aston Court  
 Kingsmead Business Park  
 Frederick Place  
 High Wycombe  
 Hp11 1JU  
 United Kingdom  
 T. +44 (0)1494 616 082  
 F. +44 (0)1494 616 243



REV	DATE	DESCRIPTION
00	11/06/2018	Emission
01	26/02/2021	Changed mould and tubes configurations