


## TPN multilayer Bag 3000ml 6-way gravity fill set with in-line filter

**CODE:** MB-1431-6W

**PCS/BOX:** 30

**WEIGHT:** 266g

1. Renolit Solmed Medipak UVO 9002 400 micron
2. Internal handle (PS)
3. Black graduated scale
4. Coex Ø4.8x6.8mm tubes (PUR-EVA)
5. Coex Ø6X8mm tube (PUR-EVA)
6. Infusion port (ABS)
7. Safety cap for infusion port (LDPE)
8. Inviolable clamp (SBC)
9. Spike port (PVC DEHP free)
10. Male connector (ABS)
11. Female connector (ABS)
12. Clamp on/off (PP)
13. PVC DEHP free Ø4.8x6.8mm 100mm tubes
14. Soft drip chamber (PVC DEHP free)
15. Tubular filter 40micron (Nylon)
16. Cover for drip chamber (PVC DEHP free)
17. 3 way Y adaptor (PVC DEHP free)
18. 4 way Y adaptor (PVC DEHP free)
19. PVC DEHP free Ø4.8x6.8mm 600mm tubes
20. 2 way spike with filter and vented capsule (ABS)
21. Female cap blister packaging (ABS)
22. Paper/PE envelope



**DEHP FREE**



### Intended use

Container for nutrition solutions prior and during administration to a patient.

### Controls

Test to check sterilization and to detect pyrogenic substances, control to check flow and seal.

Chemical control of components in accordance with the European Pharmacopoeia (Current edition).

### Sterilization

Product sterilized with ethylene oxide (E.O.). Validity: 5 years from sterilization date.

### Packing

Single package product in a Paper/PE envelope.

The product is supplied sterile and non pyrogenic.

### Storage

Normal storage procedures.

### Disposal

Please follow the current laws on hospital waste disposal.

**Miktell Ltd**  
Aston Court  
Kingsmead Business Park  
Frederick Place  
High Wycombe  
Hp11 1JU  
United Kingdom  
T. +44 (0)1494 616 082  
F. +44 (0)1494 616 243



00	29/06/2021	Emission
<b>REV</b>	<b>DATE</b>	<b>DESCRIPTION</b>
00	29/06/2021	Emission