


TPN Multilayer Bag with 6-way gravity fill set



Cod: UK0026BP - 1000 ML 30 Pcs/Box



DEHP FREE

CLASS IIA

1. TPN MULTILAYER BAG (EVA MULTILAYER TUBULAR 0,40mm WITH INTERNAL RIBS)
2. BLACK HOT PRINT
3. INFUSION PORT (ABS)
4. INVIOLENT CLAMP (SBC)
5. SPIKE PORT (PVC DEHP free)
6. MALE CONNECTOR (ABS)
7. FEMALE CONNECTOR (ABS)
8. PVC DEHP free TUBE Ø4,8x6,8mm - 30cm

9. CLAMP ON/OFF (PP)
10. 3-WAY Y ADAPTOR (PVC DEHP free)
11. PVC DEHP free TUBE Ø4,8x6,8mm - 10cm
12. 4-WAY Y ADAPTOR (PVC DEHP free)
13. PVC DEHP free TUBE Ø4,8x6,8mm - 60cm
14. 2-WAY SPIKE WITH FILTER AND VENTED CAPSULE (ABS)
15. FEMALE CAP BLISTER PACKAGING (ABS)
16. PAPER/PE ENVELOPE

INTENDED USE

Container for nutrition solutions prior and during administration to a patient. High protective barrier O₂, CO₂ and UV.

PACKING

Single package product in a Paper/PE envelope.
The product is supplied sterile and non pyrogenic.

CONTROLS

Test to check sterilization and to detect pyrogenic substances, control to check flow and seal.
Chemical control of components in accordance with the European Pharmacopoeia (Current edition).

STERILIZATION

Product sterilized with ethylene oxide (E.O.). Validity: 5 years from sterilization date.

STORAGE

Normal storage procedures.

DISPOSAL

Please follow the current laws on hospital waste disposal.

CE 0051

REV	DATE	DESCRIPTION	ISSUED BY	APPROVED BY	QA APPROVED
00	29/05/2020	Emission	L. Rossi	M. Mancarti	G. Galli
01	09/02/2021	Changed mould	L. Rossi	M. Mancarti	G. Galli