


TPN EVA Bag 5000ml with female connector

CODE: MB-1350

PCS/BOX: 40

WEIGHT: 115g

1. Eva tubular (0.35mm thickness)
2. Internal handle (PS)
3. Black graduated scale
4. Coex Ø4.8x6.8mm tubes (PUR-EVA)
5. Coex Ø6X8mm tube (PUR-EVA)
6. Infusion port (ABS)
7. Safety cap for infusion port (LDPE)
8. Inviolable clamp (SBC)
9. Spike port (PVC DEHP free)
10. Female connector (ABS)
11. Male vented cap (ABS)
12. Male cap blister packaging (ABS)
13. Paper/PE envelope



DEHP FREE



Intended use

Container for nutrition solutions prior and during administration to a patient.

Controls

Test to check sterilization and to detect pyrogenic substances, control to check flow and seal.

Chemical control of components in accordance with the European Pharmacopoeia (Current edition).

Sterilization

Product sterilized with ethylene oxide (E.O.). Validity: 5 years from sterilization date.

Packing

Single package product in a Paper/PE envelope.

The product is supplied sterile and non pyrogenic.

Storage

Normal storage procedures.

Disposal

Please follow the current laws on hospital waste disposal.

Miktell Ltd
 Aston Court
 Kingsmead Business Park
 Frederick Place
 High Wycombe
 Hp11 1JU
 United Kingdom
 T. +44 (0)1494 616 082
 F. +44 (0)1494 616 243



REV	DATE	DESCRIPTION
00	11/06/2018	Emission
01	07/06/2021	Changed mould and tubes configurations